




Digital Potentiometer 318



318
TITRATION



SYSTRONICS (INDIA) LIMITED (SYSTRONICS DIVISION)

Works : 89-92, Industrial Area, Naroda, Ahmedabad - 382 330. Gujarat, INDIA.
Phones : (079) 22813017 / 22813117, Fax : (079) 22821592,

Regd. Office : B/116-129, Supath-II Complex,
Near Juna Wadaj Bus Terminus, Ashram Road, Ahmedabad-380 013.
Ph. : 079-27557072, 27556077, Fax : 079-27552902

Email : sales@systronicsindia.com • mktg@systronicsindia.com
Website : www.systronicsindia.com

MANUFACTURERS OF ELECTRONIC INSTRUMENTS AND SCIENTIFIC EQUIPMENTS.

BRANCH ADDRESS :

BRANCHES AT : BENGALURU, BHOPAL, BHUBANESHWAR,
CHANDIGARH, CHENNAI, GUWAHATI, HYDERABAD, JAIPUR, KOLKATA,
LUCKNOW, MUMBAI, NEW DELHI, PATNA & THIRUVANANTHAPURAM.

Digital Potentiometer 318

SPECIFICATIONS

Input Configuration	Bipolar
mV	
Range	0 to ± 1999 V
Resolution	1 mV
Control	
Set Zero Range	0 to ± 1500 mV input
Stirrer Speed	0 to full speed
Switch	Hold reading when pressed
Polarization Current	10 μ A for Metal to Metal electrode
Input Impedance	10^{11} Ohms
Polarity	Automatic
Over Range	Last 3 digits blinking
Magnetic Stirrer	In-built with variable control
Display	3½ Digit LED
Power	230 V AC $\pm 10\%$, 50Hz
Dimensions	235(W) X 185(D) X 85(H) mm
Weight	1.25 Kg (Approx)
Accessories	(1) Combined Pt. & Ref. Electrode or Pt. & Ref. Electrode (2) Electrode Stand & Clamp (3) Stirrer Pedal

The unit is Injection Molded Cabinet.



*Technical Specifications, Appearance & model number are subject to change without prior notice in keeping up with the state of art.